



Long-read sequence analysis

Applications

Why long reads?




Applications

- (Genome) assembly
- Variant analysis
- Transcriptome analysis
- Epigenetics
- Metagenomics

Assembly

- Reconstructing a DNA sequence out of fragmented sequences:
 - Better understand variations between and within species
 - Reduces resources and increases accuracy for many applications
- Longer reads:
 - Less assembly errors
 - Higher contiguity (longer contigs)
 - Lower computational resources required

Assembly





LR assembly software

- Microbial:
 - HGAP
 - Trycycler
- Large(r) genomes:
 - Shasta (ONT only)
 - Raven
 - wtdbg2
 - Falcon (PacBio only)
 - Canu
 - Flye


Variant analysis

- Low accuracy long reads: not particularly good at variant analysis
- Major applications
 - Structural variation/large repeats
 - Phasing of variants



Phasing


- Variant calling:
 - DeepVariant (Pacbio + ONT)
 - Medaka (ONT)
- Phasing: Whatsp



Repeat expansion - project 2



Transcriptome analysis





Transcriptome analysis


- Full transcripts: less ambiguous identification + quantification of transcript isoforms
- Better gene model prediction (e.g. for genome annotation)
- ONT direct RNA-seq:
 - No PCR amplification
 - No reverse transcription
 - Base modification identification

Epigenetics


- Base modification in non-amplified libraries (Pacbio and ONT)
- No multiplexing
- Also possible in direct RNA-seq (ONT)




5-mC




5-hmC




5-fC




5-caC




4-mC



6-mA



8-oxoG




8-oxoA

Metagenomics

- Full length 16S genes
- Metagenome Assembled Genomes (MAGs)

New possibilities

- Cas9 based targeting (both PacBio & ONT): also base modifications
- Adaptive sampling (ONT)



<https://genomebiology.biomedcentral.com/articles/10.1186/s13059-021-02582-x>
<https://www.nature.com/articles/s42003-021-01648-7>

image: <https://store.nanoporetech.com/catalog/product/view/id/461/s/cas9-sequencing-kit/category/28/>